

CARE AND CLEANING GUIDE

Muse

Spot cleaning is always recommended. Prompt cleaning of stains will produce best results.  
Recommended cleaning steps are listed by type of stain.

REGULAR CLEANING AND MAINTENANCE

EXAMPLE

Clean the soiled area with mild soap and water, then rinse with fresh water and wipe down with a clean cloth.

Dirt  
Dust  
Grime

FOOD STAINS / OILS

Wipe affected area with a soft cloth or soft bristle brush using a none abrasive cleaner, for example Formula 409 or Fantastik. Rinse with fresh water and wipe dry.

Tea  
Ketchup  
Chocolate  
Coffee  
Juice

DIFFICULT STAINS

Clean the soiled area with a 50% solution of Isopropyl Alcohol and water. It is important to rinse the cleaned area with fresh water after applying the alcohol solution.

Eyeshadow / Mascara  
Lipstick  
Crayon  
Mustard  
Grease

HEALTHCARE / DISINFECTION

Clean the soiled area with a 10% solution of Household Bleach and water. It is important to rinse the cleaned area with fresh water after applying the bleach solution.

Blood  
Urine  
Betadine

NOTE: The information in this cleaning guide refers to performance in specific tests conducted under laboratory conditions. This information is not a guarantee and does not relieve the user from the responsibility of the proper and safe use of the product and referenced cleaning agents.